



REV	DESCRIPTION	DATE	APP'D
Α	FIRST ISSUE	10/12/12	
Α	FIRST ISSUE	10/12/12	

REVISIONS










**SECTION A-A** 



This drawing and associated designs are the property of Biopure Technology Ltd. and may not be copied or used for any purpose other than that for which they were supplied without the express written permission of:

BioPure Technology Ltd Unit M1, Hazelton Interchange Lakesmere Road Horndean Hampshire PO8 9JU

Tel: +44 (0) 23 9249 9000 Fax: +44 (0) 23 9249 9001

Email: sales@biopuretech.com Web: www.biopuretech.com

UNLESS OTHERWISE SPECIFIED: DIMENSIONS IN MM TOLERANCES: X: ±0.5 X.X: ±0.25 X.XX: ±0.1 ANGULAR: ±1°

Smooth

TITLE:

1"-1 1/2" FLAT BIO END CAP

DWG NO: PEC-100-150-FC

NAME SIGNATURE LT DRAWN CHK'D JMc

DATE 10/12/12 10/12/12 WEIGHT:

MATERIAL: Bormed HD810MO PART NO:

SCALE: NTS

PEC-100-150-FC

A4

SHEET 1 OF 1

© COPYRIGHT